

*super***Ketch**[™]

*speed***Ketch**[™]



First & original vascular
access system designed
for **single-handed** use

- A unique push-pull design
- Clean, safe and fast
- Essential for Drug Eluting Stent
- Full range, from connector to inflation kit



min*V***ASYS**

■ A unique design

- For easy, fast and accurate single handed intra vascular device insertion

■ Clean and Safe

- Combination of a haemostatic valve integrated push-pull button valve opener
- Cleaner and safer procedure with a minimized backflow thus limiting blood loss

■ Essential for DES

- Perfect insertion device respecting the integrity of the coated surface

■ Technical data

- Speed Ketch & Super Ketch 7F inner lumen compatibility
- 8 atm. pressure resistance after device introduction*

■ Ordering information

	Speed Ketch	Super Ketch	Super Ketch LSA	Torquer	Insertion Tool	Inflation device
KESU0001		20 cm side arm			Up to 0.025"	
KESU0002		20 cm side arm		0.010" to 0.035"	Up to 0.025"	
KESU0003			50 cm side arm		Up to 0.025"	
KESU0004	No side arm				Up to 0.025"	
KESU0005			50 cm side arm	0.010" to 0.035"	Up to 0.025"	
KESU0006	No side arm			0.010" to 0.035"	Up to 0.025"	
KISK0001		20 cm side arm		0.010" to 0.035"	Up to 0.025"	30 atm / 20 cc

(*) Haemostatic valve resistance measured in vitro, by applying injection pressure through the side port.

